


<b>MODEL</b>	<b>MIO-VEIN</b>	
<b>VERSION</b>	<b>FES ELECTROTHERAPY</b>	
<p style="text-align: center;"><b>CONTENTS:</b></p> <ul style="list-style-type: none"> <li>• MIO-VEIN</li> <li>• 2 outputs</li> <li>• 2 connection cables</li> <li>• 4 oval electrodes for calf</li> <li>• Belt clip</li> <li>• Ni-Mh internal rechargeable battery with LL (long-lasting) technology.</li> <li>• Battery charger</li> <li>• User manual</li> <li>• Carriage bag</li> </ul>		
<p><b>MAIN TECHNICAL FEATURES</b></p> <ul style="list-style-type: none"> <li>• Two independent channels for the connection of 4 electrodes;</li> <li>• easy to use digital keyboard;</li> <li>• adjustable intensity</li> <li>• on-line guide on display for fast instructions;</li> <li>• reflective LCD display;</li> <li>• FES specific program for venous stimulation;</li> <li>• operating with internal rechargeable battery;</li> <li>• 60 mApp maximum intensity;</li> <li>• width impulse 30-240 <math>\mu</math>s;</li> <li>• CE0476 medical device</li> </ul>		
<p><b>WHY A CALF STIMULATOR?</b></p> <p><b>Mio-Vein</b> is the new 2 channels FES device (Functional Electrical Stimulation) able to stimulate calf muscles and helps the correct venous circulation.</p> <p>Its stimulation reactivates calf muscular pump that is mainly responsible of venous return. Calf muscles contraction promotes the compression of superficial and deep veins. As a consequence blood reaches easily the heart.</p> <p><b>Mio-Vein</b> can be used not only to <b>help the natural blood</b> flow but also to <b>prevent venous insufficiency</b> and <b>deep vein thrombosis</b>. Moreover it can be used in presence of lower limb oedemas, venous ulcers and venous diseases.</p> <p><b>Mio-Vein</b> represents an important support to compressive stocking therapy and finally to venous stasis treatment.</p>		
<p><b>Size:</b> 114 x 68 x 26 mm</p>		